

PH522 Intro: Gaussian Beams

Objective: You will measure the variation with distance of a Gaussian beam, and compare with theory. You will also study the transformation of a Gaussian beam by a lens and check with ABCD matrices.

Equipment:

- 1) Commercial 2 mW HeNe laser
- 2) a positive lens of focal length 30 cm.
- 2) Power meter
- 3) Knife edge on a translation stage

Procedure:

- 1) Observe the HeNe laser spot on the wall and note (without measurement) that it is much larger than at the laser output.
- 2) Position the power meter far downstream so that it intercepts the full laser power.
- 3) Position the translation stage and knife edge near the laser output, and record the measured power as a function of the position of the knife edge. Use your data to determine the beam waist, $w(z)$.
- 4) Repeat this procedure for at least three other positions z . Make sure you include some measurements very far away from the laser output.
- 5) Plot your measured values of $w(z)$ vs. the positions z at which they were measured. You might also choose to plot $w(z)^2$ vs. z^2 , which should be linear.
- 6) Assume that w_0 is at the HeNe output and check if $w(z)$ follows the formula for Gaussian beam propagation.
- 7) Place a 30 cm lens 1 m in front of the laser. After the lens, the beam will have a new focal position and size at that position.
- 8) Using an index card, try to estimate the location of the beam waist by observing the variation of the spot size. Record your best estimate of the new w_0 location in your lab book. Compare with predictions from the ABCD rule.
- 9) Measure the asymptotic angular divergence of the new beam by observing the spot a large distance beyond the lens and assuming that the diameter of the spot is $2w(z)$. With this, estimate the new value of w_0 .